pH Buffer Solutions

• Safety Data Sheets

 Safety data sheets for all Hanna solutions are available at hannainst.com or upon request.

Expiration date

 The production batch number, expiration date, and temperature correlation table are reported on all Hanna calibration solutions.

NIST traceability

 Standardized using a pH meter calibrated by means of two standard solutions prepared from NIST standard reference materials.



Air-tight bottles

 Air-tight bottle with tamper-proof seal of freshness to ensure quality.

• Single use sachets

 Opaque packaging prevents oxidation from UV light that could alter the value. Every sachet is as fresh as the day it was packaged.

• FDA compliant bottles (HI80xx)

 Hanna solutions are offered in opaque bottles that meet FDA requirements.

4.01 pH Buffer Solution

This buffer value is widely used in water purification plants, in the food industry, and wherever the pH is expected to be slightly acidic.





4.01 pH @ 25°C - Bottles

Code	Size	FDA Bottle	Certificate of Analysis
HI7004/1G	1 G (3.78 L) (color coded solution)	on request	
HI7004/1L	1 L (color coded solution)	on request	
HI7004L	500 mL	on request	
HI7004L/C	500 mL	•	
HI7004C	500 mL (color coded solution)	on request	
HI7004M	230 mL	on request	
HI7004-050	500 mL (GroLine®)	•	
HI7004-023	230 mL (GroLine)	•	
HI7004-012	120 mL (GroLine)	•	
HI8004L	500 mL	•	•
HI8004L/C	500 mL	•	•

4.01 pH @ 25°C - Sachets

Code	Size	Package	Certificate of Analysis
HI70004C	20 mL	25 pcs.	•
HI70004G	20 mL (GroLine)	25 pcs.	•
HI70004P	20 mL	25 pcs.	
HI700044P	20 mL (Pool Line)	25 pcs.	

4.01 and 7.01 pH @ 25°C - Sachets

Code	Size	Package	Certificate of Analysis
HI77400C	20 mL	10 pcs., 5 ea	•
HI77400P	20 mL	10 pcs., 5 ea	

